



Anti-IER2 monoclonal antibody (DCABH-695)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Product Overview | Mouse monoclonal to IER2 |
| Antigen Description | Mouse monoclonal antibody raised against a full-length recombinant IER2. |
| Immunogen | Recombinant full length protein, corresponding to amino acids 1-224 of Human IER2 with a proprietary tag (AAH03625). |
| Isotype | IgG2a |
| Source/Host | Mouse |
| Species Reactivity | Human |
| Conjugate | Unconjugated |
| Applications | WB |
| Format | Liquid |
| Size | 100 µg |
| Buffer | pH: 7.20; Constituent: 99% PBS |
| Preservative | None |
| Storage | Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing. |

GENE INFORMATION

| | |
|------------------------|--|
| Gene Name | IER2 immediate early response 2 [Homo sapiens] |
| Official Symbol | IER2 |

| | |
|----------------------------|---|
| Synonyms | IER2; immediate early response 2; immediate early response gene 2 protein; ETR101; MGC15265; MGC111472; |
| Entrez Gene ID | 9592 |
| Protein Refseq | NP_004898 |
| UniProt ID | A0A024R7H1 |
| Chromosome Location | 19p13.13 |