



Anti-SPSB1 monoclonal antibody (DCABH-702)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to SPSB1
Antigen Description	Probable substrate recognition component of a SCF-like ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.
Immunogen	Recombinant full length protein, corresponding to amino acids 1-274 of Human SPSB1 with a proprietary tag (AAH15711).
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA, Sandwich ELISA
Positive Control	Recombinant Protein.
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Ship	Shipped at 4°C.

GENE INFORMATION

Gene Name	SPSB1 splA/ryanodine receptor domain and SOCS box containing 1 [Homo sapiens]
Official Symbol	SPSB1
Synonyms	SPSB1; splA/ryanodine receptor domain and SOCS box containing 1; SPRY domain-containing SOCS box protein 1; SSB 1; novel SPRY domain containing protein; SPRY domain-containing SOCS box protein SSB-1; SSB1; SSB-1;
Entrez Gene ID	80176
Protein Refseq	NP_079382
UniProt ID	A0A024R4G8
Chromosome Location	1p36.22
Pathway	Adaptive Immune System, organism-specific biosystem; Antigen processing: Ubiquitination & Proteasome degradation, organism-specific biosystem; Class I MHC mediated antigen processing & presentation, organism-specific biosystem; Immune System, organism-specific biosystem;