



# Anti-FAM127C monoclonal antibody, clone 3E5 (DCABH-746)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse monoclonal to FAM127C
<b>Immunogen</b>	Recombinant full length protein, corresponding to amino acids 1-113 of Human FAM127C produced in HEK293T cells (NP_001071641).
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	3E5
<b>Purification</b>	Purified from Mouse ascites fluids by affinity chromatography.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF
<b>Positive Control</b>	FAM127C transfected HEK293T cell lysate and COS7 cells transiently transfected with FAM127C
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.30; Preservative: 0.02% Sodium azide; Constituents: 48% PBS, 1% BSA, 50% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	store at -20°C. Avoid freeze / thaw cycles.

Ship

Shipped at 4°C.

---

## GENE INFORMATION

Gene Name	<a href="#">FAM127C family with sequence similarity 127, member C [ Homo sapiens ]</a>
Official Symbol	FAM127C
Synonyms	FAM127C; family with sequence similarity 127, member C; protein FAM127C; CXX1c; MAR8B; mammalian retrotransposon derived protein 8B; FLJ25577;
Entrez Gene ID	<a href="#">441518</a>
Protein Refseq	<a href="#">NP_001071641</a>
UniProt ID	<a href="#">Q17RB0</a>
Chromosome Location	Xq26.3