



Anti-NECAB1 monoclonal antibody, clone 5H9 (DCABH-779)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to NECAB1
Immunogen	Recombinant full length Human NECAB1 (NP_071746) produced in HEK293T cells.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5H9
Purification	Purified from mouse ascites fluids by affinity chromatography.
Conjugate	Unconjugated
Applications	WB, Flow Cyt
Positive Control	HEK293T cells transfected with pCMV6-ENTRY NECAB1 cDNA, Hela and Jurkat cells
Format	Liquid
Size	100 µl
Buffer	pH: 7.30; Preservative: 0.02% Sodium azide; Constituents: 48% PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	store at -20°C. Avoid repeated freeze / thaw cycles.
Ship	Shipped at 4°C.

GENE INFORMATION

Gene Name	NECAB1 N-terminal EF-hand calcium binding protein 1 [Homo sapiens]
Official Symbol	NECAB1
Synonyms	NECAB1; N-terminal EF-hand calcium binding protein 1; EF hand calcium binding protein 1 , EFCBP1; N-terminal EF-hand calcium-binding protein 1; neuronal calcium binding protein; EF hand calcium binding protein 1; EF-hand calcium binding protein 1; EF-hand
Entrez Gene ID	64168
Protein Refseq	NP_071746
UniProt ID	Q8N987
Chromosome Location	8q21.3
Function	calcium ion binding;