



Anti-CCNB1IP1 monoclonal antibody, clone 5I4 (DCABH-844)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to CCNB1IP1
Antigen Description	E3 ubiquitin-protein ligase. Modulates cyclin B levels and participates in the regulation of cell cycle progression through the G2 phase. Overexpression causes delayed entry into mitosis.
Immunogen	Recombinant full length Human CCNB1IP1 produced in E.coli (NP_878269).
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	5I4
Purification	This antibody is purified from Mouse ascites fluid by affinity chromatography.
Conjugate	Unconjugated
Applications	WB, Flow Cyt
Positive Control	HEK293T cell lysate transfected with pCMV6-ENTRY CCNB1IP1 cDNA; HeLa and Jurkat cells.
Format	Liquid
Size	100 µl
Buffer	pH: 7.30; Preservative: 0.02% Sodium azide; Constituents: 48% PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide

Storage	store at -20°C. Avoid repeated freeze / thaw cycles.
----------------	--

Ship	Shipped at 4°C.
-------------	-----------------

GENE INFORMATION

Gene Name	CCNB1IP1 cyclin B1 interacting protein 1, E3 ubiquitin protein ligase [Homo sapiens]
------------------	--

Official Symbol	CCNB1IP1
------------------------	----------

Synonyms	CCNB1IP1; cyclin B1 interacting protein 1, E3 ubiquitin protein ligase; C14orf18, chromosome 14 open reading frame 18 , cyclin B1 interacting protein 1; E3 ubiquitin-protein ligase CCNB1IP1; HEI10; human enhancer of invasion 10; enhancer of invasion 10;
-----------------	--

Entrez Gene ID	57820
-----------------------	-----------------------

Protein Refseq	NP_067001
-----------------------	---------------------------

UniProt ID	Q9NPC3
-------------------	------------------------

Chromosome Location	14q11.2
----------------------------	---------

Function	ligase activity; metal ion binding; protein binding;
-----------------	--