



Anti-HBA1 monoclonal antibody, clone 8F2G (DCABH-9765)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to Hemoglobin
Antigen Description	Involved in oxygen transport from the lung to the various peripheral tissues.
Immunogen	Human plasma
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	8F2G
Conjugate	Unconjugated
Applications	ELISA, WB, Flow Cyt
Positive Control	Human plasma.
Format	Liquid
Size	100 µl
Buffer	Preservative: 0.03% Sodium Azide; Constituents: 50% Glycerol, 0.01% BSA, HEPES, 0.15M Sodium chloride
Preservative	0.03% Sodium Azide
Storage	store at -20°C. Avoid freeze / thaw cycles.

Ship

Shipped at 4°C.

GENE INFORMATION

Gene Name	HBA1 hemoglobin, alpha 1 [Homo sapiens]
Official Symbol	HBA1
Synonyms	HBA1; hemoglobin, alpha 1; hemoglobin subunit alpha; alpha-globin; alpha-1 globin; alpha-1-globin; alpha one globin; hemoglobin alpha chain; hemoglobin alpha-1 chain; hemoglobin alpha 1 globin chain; HBH; CD31; MGC126895; MGC126897;
Entrez Gene ID	3039
Protein Refseq	NP_000549
UniProt ID	D1MGQ2
Chromosome Location	16p13.3
Pathway	African trypanosomiasis, organism-specific biosystem; African trypanosomiasis, conserved biosystem; Malaria, organism-specific biosystem; Malaria, conserved biosystem; Selenium Pathway, organism-specific biosystem;