



Anti-CCL26 monoclonal antibody, clone 226003 (DCABY-4018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCL26 (chemokine (C-C motif) ligand 26) is a protein-coding gene. Diseases associated with CCL26 include eosinophilic esophagitis, and heart tumor. GO annotations related to this gene include chemokine activity. An important paralog of this gene is CCL15.
Specificity	Detects human CCL26/Eotaxin-3 in ELISAs and Western blots. Does not cross-react with recombinant human (rh) CCL11/Eotaxin, recombinant mouse (rm) CCL11/Eotaxin, rhCCL24/Eotaxin-2, rmCCL24/Eotaxin-2, rhCCL3, rhCCL4, rmCCL3, rmCCL4, rhCCL5, rhCCL13, or rhCCL14.
Immunogen	E. coli-derived recombinant human CCL26/Eotaxin-3. Thr24-Leu94 Accession Number Q9Y258
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	226003
Purification	Protein A or G purified from ascites
Conjugate	Unconjugated
Applications	Western Blot, ELISA Capture (Matched Pair)
Format	Liquid
Size	500 µg
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose.

Preservative	None
Storage	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <p>12 months from date of receipt, -20 to -70 °C as supplied.</p> <p>1 month from date of receipt, 2 to 8 °C, reconstituted.</p> <p>6 months from date of receipt, -20 to -70 °C, reconstituted.</p>

GENE INFORMATION

Gene Name	CCL26 chemokine (C-C motif) ligand 26 [Homo sapiens (human)]
Official Symbol	CCL26
Synonyms	CCL26; chemokine (C-C motif) ligand 26; IMAC; TSC-1; MIP-4a; SCYA26; MIP-4alpha; C-C motif chemokine 26; eotaxin-3; MIP-4-alpha; chemokine N1; CC chemokine IMAC; thymic stroma chemokine-1; small inducible cytokine A26; small-inducible cytokine A26; macrop
Entrez Gene ID	10344
Protein Refseq	NP_006063
UniProt ID	Q9Y258
Chromosome Location	7q11.23
Pathway	Chemokine signaling pathway; Cytokine-cytokine receptor interaction; IL4-mediated signaling events;
Function	chemokine activity;