



Magic™ Anti-Parvovirus Monoclonal antibody, Clone N29723QW2 (DCABY-4590)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Dog
Target	Parvovirus
Immunogen	Parvovirus antibody was raised in mouse using canine parvovirus as the immunogen.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Parvovirus
Clone	N29723QW2
Purification	Parvovirus antibody was purified by chromatography on protein A Sepharose.
Conjugate	Unconjugated
Applications	ELISA Pr*, Immunoassays, Immunodiffusion, Lateral Flow Pr* Suggested pair for sandwich ELISA (Capture - Detection): DCABY-4590 - DCABY-4591
Format	Liquid
Size	1 mg
Buffer	IgG purified by protein A chromatography and supplied in PBS, with 0.05% NaN ₃ .
Preservative	0.05% Sodium Azide
Storage	Store at 4 °C

BACKGROUND

Introduction	Parvoviruses are typically linear, non-segmented single-stranded DNA viruses, with an average genome size of 5000 nucleotides. Parvoviruses are some of the smallest viruses (hence the name, from Latin parvus meaning small) and are 18–26 nm in diameter.
Keywords	Parvovirus; Canine parvovirus;