



## User's Manual

# C-reactive Protein ELISA Kit

**REF** DEIA012J

 96T

**RUO**

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The assay is a sandwich enzyme immuno assay intended for the quantitative determination of C-reactive Protein in dried blood samples.

### Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

### Precision

Intra-Assay Precision: 5.5%

Inter-Assay Precision: 11.6%-13.8%

### Sensitivity

0.921 ng/mL