



User's Manual

PMN-Elastase (1-point-calibr.) ELISA Kit



DEIA033J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The described Enzyme-Linked-Immuno-Sorbent-assay (ELISA) is intended for the quantitative determination of PMN Elastase in serum, and seminal plasma.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Precision

Serum

Intra-Assay Precision: 0.13%-0.5%

Inter-Assay Precision: 5.45%-6.86%

Plasma

Intra-Assay Precision: 1.98%-2.68%

Inter-Assay Precision: 6.19-11.3%

Sensitivity

0.011 ng/mL