



## User's Manual

# Calprotectin (MRP 8/14) (1-point-calibration) ELISA Kit

REF

DEIA044J



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)



Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The described Calprotectin ELISA is intended for the quantitative determination of calprotectin (MRP (8/14, S100A8/A9) in urine.

### Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

### Precision

Intra-Assay Precision: 3.2%-5.6%

Inter-Assay Precision: 4.4%-8.9%

### Sensitivity

2.099 ng/mL

### Specificity

No cross reactivity was observed to the following plasma proteins: Lysozyme; PMN-Elastase; Myeloperoxidase; Lactoferrin.