



User's Manual

Thiamphenicol drug residue ELISA Kit

REF

DEIANJ33



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Thiamphenicol drug residue ELISA Kit is intended for the quantitative determination of Thiamphenicol drug residue in raw milk, egg, seafood, tissue.

Reagents And Materials Provided

1. Microtiter well plate: 8 microwells/strip, 12 strips/plate.
2. Standards
3. Washing Solution
4. Enzyme Conjugate Diluent
5. Substrate Solution (TMB)
6. Stop Solution

Storage

Store all reagents at 2-8°C. Do not freeze. Reagents may be used until the expiration date on the box.