



Trimethoprim Standard (98%) (DISNJ03)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Expression System	N/A
Purity	98%
Conjugate	N/A
Cas.No	738-70-5
Applications	ELISA
Molecular Weight	290.32 g/mol
Format	Solid
Size	5 g
Preservative	None
Storage	Room temperature

BACKGROUND

Introduction	Trimethoprim (TMP) is an antibiotic used mainly in the treatment of bladder infections.
Keywords	Trimethoprim; Antibacterial; Antiprotozoal