



## Doxycycline Standard (97%) (DISNJ25)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Purity	97%
Conjugate	N/A
Cas.No	564-25-0
Applications	ELISA
Molecular Weight	444.43 g/mol
Format	Solid
Size	25 mg
Preservative	None
Storage	Room temperature

### BACKGROUND

Introduction	Doxycycline is an antibiotic that is used in the treatment of a number of types of infections caused by bacteria and protozoa. It is useful for bacterial pneumonia, acne, chlamydia infections, early Lyme disease, cholera and syphilis.
Keywords	Doxycycline; Tetracycline; Antibiotic