



Magic™ Anti-IFN-γ Monoclonal antibody, Clone 2B7.2G7 (DMAB-S004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	cross reacts with cat, cougar
Immunogen	recombinant lion IFN-γ (expressed in HEK293 cells)
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	lion
Purification	GammaBindPlus purified
Conjugate	Unconjugated
Applications	ELISA(Cap), WB We recommend the following for sandwich ELISA (Capture - Detection): DMAB-S001 - DMAB-S002, DMAB-S003
Format	Liquid
Concentration	0.6 mg/ml
Size	100 µg
Buffer	PBS
Preservative	None
Storage	Store at 4° C short term (1-2 weeks) Aliquot and store at -21° C