



Anti-Methotrexate monoclonal antibody (DMAB-WZ0004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	N/A
Immunogen	MTX-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MTX04
Conjugate	Unconjugated
Applications	N/A
Size	500 µg
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	METHOTREXATE; METHYLAMINOPTERIN; METHYLAMINOPTERINE; ANTIFOLAN; L-(+)-AMETHOPTERIN; 4-AMINO-N 10-METHYLPTEROYL-L-GLUTAMIC ACID
----------	--