



Mouse Anti-Benzodiazepine monoclonal antibody, clone H12392N (DMAB3039)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	MAb to Benzodiazepine. Monoclonal Antibody to Benzodiazepine
Specificity	Benzodiazepine
Target	Benzodiazepine
Immunogen	Benzodiazepine [BSA].
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	H12392N
Purification	> 95% pure (SDS-PAGE). Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.02% Sodium Azide

Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles
----------------	--

BACKGROUND

Keywords	BZD
-----------------	-----

GENE INFORMATION

Gene Name	BZD
------------------	-----

Synonyms	BZD
-----------------	-----
