



Mouse Anti-Buprenorphine monoclonal antibody, clone C2702N (DMAB3069)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	MAb to Buprenorphine. Monoclonal Antibody to Buprenorphine
Specificity	Buprenorphine
Target	Buprenorphine
Immunogen	Buprenorphine [BSA].
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	C2702N
Purification	> 95% pure (SDS-PAGE). Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM PBS, 50 mM Sodium Chloride, pH 8.0.
Preservative	0.1% Sodium Azide

Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles
----------------	--

BACKGROUND

Keywords	AC1LD8IH; Buprenex; Buprenorfina; Buprenorphinum; DEA No. 9064; Norbuprenorphine; Subutex; Temgesic;
-----------------	--

GENE INFORMATION

Gene Name	AC1LD8IH
Synonyms	AC1LD8IH; Buprenex; Buprenorfina; Buprenorphinum; DEA No. 9064; Norbuprenorphine; Subutex; Temgesic;
