



Mouse Anti-Methadone monoclonal antibody, clone C2791N (DMAB3910)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	MAb to Methadone. Monoclonal Antibody to Methadone
Specificity	Methadone
Target	Methadone
Immunogen	Methadone [BSA].
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	C2791N
Purification	> 95% pure (SDS-PAGE). Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.02% Sodium Azide

Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles
----------------	--

BACKGROUND

Keywords	Methadone;
-----------------	------------

GENE INFORMATION

Gene Name	Methadone
Synonyms	Methadone;
