



Mouse Anti-Chicken IgY Monoclonal Antibody, clone H-2 [FITC] (DMAB4621)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Chicken
Clone	H-2
Purification	Purified
Conjugate	FITC
Applications	ELISA, FLISA, FC, IHC, IP
Format	Liquid
Concentration	Lot specific
Size	0.5 mg
Buffer	Phosphate buffered saline containing 0.1% sodium azide.
Preservative	0.1% sodium azide
Storage	Store at 2-8°C; Avoid exposure to light.
Ship	Wet ice
Warnings	This product is for research use only and is not intended for diagnostic use.

BACKGROUND

Introduction	IgY is a type of immunoglobulin which is the major antibody in bird, reptile, and lungfish blood. It is also found in high concentrations in chicken egg yolk.
Keywords	IgY Fc;IgY;Immunoglobulin Y;IgY, Heavy Chain;IgY Fab