



# Magic™ Anti-CD15 monoclonal antibody, clone Darb-4 (DMAB4941)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	MAb to CD15; Mouse Monoclonal Antibody to Human CD15
<b>Specificity</b>	In Western blotting of HL60 cell lysates the antibody labels a major band of ~220 kD corresponding to the expected molecular weight of CD15(6).
<b>Immunogen</b>	Immuno-affinity purified X-Hapten oligopeptide fragment from myelomonocytic leukemia cells
<b>Isotype</b>	IgM
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	Darb-4
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC
<b>Format</b>	Monoclonal Mouse antibody provided in liquid form as cell culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 0.015 mol/L sodium azide. This product contains stabilizing protein.
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Store at 2-8°C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens.

# GENE INFORMATION

Entrez Gene ID [2526](#)

---

UniProt ID [P22083](#)

---