



# Mouse anti-Cyanazine Monoclonal antibody, clone IZC 270-02 (DMAB5751)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Mouse monoclonal to cyanazine.
Antigen Description	The pesticide (herbicide) cyanazine.
Isotype	IgG2a
Source/Host	Mouse
Clone	IZC 270-02
Conjugate	Unconjugated
Applications	ELISA
Format	0.01 M PB, pH 7.4, with 0.5 M NaCl.
Concentration	1 mg/mL, based on E(1%)=14.0 at A280
Storage	2-8°C.

## BACKGROUND

Introduction	Cyanazine is a triazine herbicide used to control grasses and broadleaf weeds. It is now classified as restricted use as it is thought to be teratogenic.
Keywords	Cyanazine; cyanazine; Bladex; DW3418; Fortrol; Payze; 2-(4-chloro-6-ethylamino-1,3,5-triazin-2-ylamino)-2-methylpropionitrile