



Mouse Anti-Morphine monoclonal antibody, clone C2811N (DMAB7248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	MAb to Morphine. Monoclonal Antibody to Morphine
Specificity	Morphine
Target	Morphine
Immunogen	Morphine [BSA]
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	C2811N
Purification	> 95% pure (SDS-PAGE). Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM PBS, 50 mM Sodium Chloride, pH 8.0.
Preservative	0.1% Sodium Azide

Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles
---------	--

BACKGROUND

Keywords	Morphine
----------	----------

GENE INFORMATION

Gene Name	Morphine
Synonyms	Morphine
