



# Goat Anti Chicken IgY Fc Polyclonal Antibody (DMAB8564)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Highly purified chicken IgY Fc fragment isolated from normal chicken serum IgY.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Chicken
<b>Purification</b>	Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB, IHC, IP Recommended concentration: ELISA: 2 ug/ml, IP: 2 ug/mg, WB: 0.5-0.05 ug/ml
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	Phosphate buffered saline, pH 7.2 with 0.1% sodium azide.
<b>Preservative</b>	0.1% sodium azide
<b>Storage</b>	Store at 4°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Introduction** IgY is a type of immunoglobulin which is the major antibody in bird, reptile, and lungfish blood. It is also found in high concentrations in chicken egg yolk.

---

**Keywords** IgY Fc;IgY;Immunoglobulin Y;IgY, Heavy Chain;IgY Fab

---