



Donkey Anti-Chicken IgY (H+L) Polyclonal Antibody (DMAB8589)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified chicken IgY.
Isotype	IgG
Source/Host	Donkey
Species Reactivity	Chicken
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	Lot specific
Size	0.5 mg
Buffer	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide.
Preservative	0.075% sodium azide
Storage	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction IgY is a type of immunoglobulin which is the major antibody in bird, reptile, and lungfish blood. It is also found in high concentrations in chicken egg yolk.

Keywords IgY;Immunoglobulin Y;IgY Fc;IgY, Heavy Chain;IgY Fab;IgY h+l
