



Magic™ Anti-Chlamydia trachomatis Monoclonal Antibody, clone IN042 (DMABP- L05)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Chlamydia LPS
Target	C. trachomatis
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	C. trachomatis
Clone	IN042
Conjugate	Unconjugated
Applications	LFIA, ELISA
Format	Protein A Purified Monoclonal Antibody
Size	100 µg, 1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Avoid repeated freezing and thawing.