



Magic™ Anti-Human Immunoglobulin E (IgE) Monoclonal Antibody, clone GB428 (DMABP-L13)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Target | IgE |
| Immunogen | Human IgE |
| Isotype | IgG2b |
| Source/Host | Mouse |
| Species Reactivity | Human |
| Clone | GB428 |
| Conjugate | Unconjugated |
| Applications | ELISA |
| Format | Purified Monoclonal Antibody |
| Size | 1 mg |
| Buffer | 0.05 M Phosphate Buffered Saline, pH 7.4 |
| Preservative | 0.1% Sodium Azide |
| Storage | 2-8°C |