



# Magic™ Anti-Human Immunoglobulin E (IgE) Monoclonal Antibody, clone GB428 (DMABP- L13)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	IgE
<b>Immunogen</b>	Human IgE
<b>Isotype</b>	IgG2b
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	GB428
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Purified Monoclonal Antibody
<b>Size</b>	1 mg
<b>Buffer</b>	0.05 M Phosphate Buffered Saline, pH 7.4
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	2-8°C