



Magic™ Anti-Listeria monocytogenes Monoclonal Antibody, clone IN329 (DMABP- L16)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	L. monocytogenes p60
Target	L. monocytogenes
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	L. monocytogenes
Clone	IN329
Conjugate	Unconjugated
Applications	ELISA, IFA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Short Term: 2-8°C Long Term: -20°C