



Magic™ Anti-Luteinizing Hormone (LH) Monoclonal Antibody, clone GB424 (DMABP-L18)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	LH
Immunogen	Human Luteinizing hormone β -subunit
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	GB424
Conjugate	Unconjugated
Applications	ELISA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.05 M Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	2-8°C
Warnings	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

GENE INFORMATION

Entrez Gene ID [403959](#)

UniProt ID [P18842](#)
