



Magic™ Anti-Malaria (HRP-2) Monoclonal Antibody, clone IN097 (DMABP-L21)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	P. falciparum LDH
Immunogen	Recombinant Malaria HRP-2 antigen
Source/Host	Mouse
Species Reactivity	P. falciparum
Clone	IN097
Conjugate	Unconjugated
Applications	LFIA, ELISA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	None
Storage	Short Term (≤ 2 weeks): Refrigerated at 2-8°C. Long Term: Frozen at -20°C. Avoid repeated freeze/thaw cycles.