



# Magic™ Anti-Myoglobin (Mb) Monoclonal Antibody, clone IN323 (DMABP-L25)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	MB
<b>Immunogen</b>	Native Human Myoglobin
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	IN323
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	LFIA
<b>Format</b>	Liquid; Purified Monoclonal Antibody
<b>Size</b>	1 mg
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing
<b>Warnings</b>	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

## GENE INFORMATION

Entrez Gene ID [4151](#)

---

UniProt ID [P02144](#)

---