



Magic™ Anti-Norovirus genogroup I Monoclonal Antibody, clone IN033 (DMABP-L26)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Group I capsid protein
Target	NOV
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	N/A
Clone	IN033
Conjugate	Unconjugated
Applications	LFIA, ELISA, WB
Format	Protein A Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Avoid repeated freezing and thawing.
Warnings	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

GENE INFORMATION

Entrez Gene ID [4856](#)

UniProt ID [A0A024R9J4](#)
