



# Magic™ Anti-NT-proBNP (aa 13-27) Monoclonal Antibody, clone IN320 (DMABP-L30)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	epitope 13-27 aa
Target	NPPB
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	IN320
Conjugate	Unconjugated
Applications	LFIA
Format	Liquid; Purified Monoclonal Antibody
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	None
Storage	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C . Avoid repeated freezing and thawing

## GENE INFORMATION

Entrez Gene ID [4879](#)

UniProt ID

[P16860](#)

---