



# Magic™ Anti-NT-proBNP (aa 38-44) Monoclonal Antibody, clone IN257 (DMABP-L34)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	38-44 aa epitope
Target	NPPB
Immunogen	Native Human ProBNP antigen
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	IN257
Conjugate	Unconjugated
Applications	LFIA
Format	Liquid; Purified Monoclonal Antibody
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	None
Storage	Short Term ( $\leq 2$ weeks): Refrigerated at 2-8°C. Long Term: Frozen at -20°C. Avoid repeated freeze/thaw cycles.

## GENE INFORMATION

Entrez Gene ID [4879](#)

---

UniProt ID [P16860](#)

---