



Magic™ Anti-Parvovirus, Canine Monoclonal Antibody, clone IN204 (DMABP-L37)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Canine parvovirus
Target	CPV
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Dog
Clone	IN204
Conjugate	Unconjugated
Applications	LFIA, ELISA
Format	Protein A Purified Monoclonal Antibody
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	2-8°C
Warnings	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.