



Magic™ Anti-Respiratory Syncytial Virus (RSV) Monoclonal Antibody, clone K470 (DMABP-L44)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Specificity | Fusion protein of RSV types A & B |
| Target | RSV |
| Isotype | IgG2a |
| Source/Host | Mouse |
| Species Reactivity | N/A |
| Clone | K470 |
| Conjugate | Unconjugated |
| Applications | LFIA, ELISA, IFA |
| Format | Protein A Purified Monoclonal Antibody |
| Size | 1 mg |
| Buffer | 0.01 M Phosphate Buffered Saline, pH 7.2 |
| Preservative | 0.1% Sodium Azide |
| Storage | Store at -20°C. Avoid repeated freezing and thawing. |
| Warnings | This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. |