



Magic™ Anti-*Streptococcus pneumoniae* Monoclonal Antibody, clone IN299 (DMABP- L47)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	<i>S. pneumoniae</i> common cell wall polysaccharide (urinary antigen)
Target	<i>S. pneumoniae</i>
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. pneumoniae</i>
Clone	IN299
Conjugate	Unconjugated
Applications	LFIA, ELISA, IFA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	2-8°C
Warnings	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.