



Magic™ Anti-Thyroid Stimulating Hormone (TSH) Monoclonal Antibody, clone GB427 (DMABP-L50)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|---|
| Target | TSH |
| Immunogen | Human TSH β -subunit |
| Isotype | IgG |
| Source/Host | Mouse |
| Species Reactivity | Human |
| Clone | GB427 |
| Conjugate | Unconjugated |
| Applications | ELISA |
| Format | Purified Monoclonal Antibody |
| Size | 1 mg |
| Buffer | 0.05 M Phosphate Buffered Saline pH 7.4 |
| Preservative | 0.1% Sodium Azide |
| Storage | Short Term (≤ 2 weeks): Refrigerated at 2-8°C. Long Term: Frozen at -20°C. Avoid repeated freeze/thaw cycles. |

GENE INFORMATION

Entrez Gene ID [1081](#)

UniProt ID [P01215](#)
