



# Magic™ Anti-Thyroid Stimulating Hormone (TSH) Monoclonal Antibody, clone IN224 [HRP] (DMABP-L51)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	TSH
Immunogen	Human TSH $\beta$ -subunit
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	IN224
Conjugate	HRP
Applications	ELISA
Format	Purified Monoclonal Antibody-HRP Conjugate
Size	1 mg
Buffer	protein stabilization buffer, pH 7.6
Preservative	None
Storage	4-8°C

## GENE INFORMATION

Entrez Gene ID [1081](#)

---

UniProt ID [P01215](#)

---