



# Anti-Astrovirus Monoclonal antibody, Clone 90000000 (DMABT-51096MA)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Anti AstroVirusMouse Anti AstroVirus
<b>Target</b>	Astrovirus
<b>Immunogen</b>	Astrovirus type 2 expressed in LLC MK2 cells.
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Astrovirus
<b>Clone</b>	90000000
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, ELISA, FC, IP, WB
<b>Format</b>	Purified IgG - liquid
<b>Concentration</b>	IgG concentration 1.0mg/ml
<b>Size</b>	1 mg
<b>Buffer</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4 °C or at -20 °C if preferred.Storage in frost-free freezers is not recommended.This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation

before use.

---

## BACKGROUND

### Introduction

A genus of small, circular RNA viruses in the family ASTROVIRIDAE. They cause gastroenteritis and are found in the stools of several vertebrates including humans. Transmission is by the fecal-oral route. There are at least eight human serotypes.

---

### Keywords

Astrovirus; Mamastrovirus; Avastrovirus

---