



Anti-FOXL1 monoclonal antibody, clone 3G7 (DMABT-H13781)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-FOXL1 Monoclonal Antibody
Target	FOXL1
Immunogen	FOXL1 (NP_005241, 132 a.a. ~ 240 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3G7
Conjugate	Unconjugated
Applications	WB, IHC, IF, ELISA
Sequence Similarities	MFENGNYRRRKRPKPGGAPEAKRPRAETHQRSAAEQPEAGSGAGGSGPAISRLQAAPA GPSPLLDGPSPPAPLHWPGTASPNEDAGDAAQGAAAVAVGQAARTGDGP
Buffer	In 1x PBS, pH 7.2
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FOXL1 forkhead box L1 [Homo sapiens]
Official Symbol	FOXL1
Synonyms	FOXL1; forkhead box L1; FKHL11; forkhead box protein L1; FKH6; FREAC7; FREAC-7; forkhead-like 11; forkhead-related protein FKHL11; forkhead-related transcription factor 7;
Entrez Gene ID	2300
Protein Refseq	NP_005241
UniProt ID	Q12952
Chromosome Location	16q24
Function	DNA binding; DNA binding, bending; DNA binding, bending; double-stranded DNA binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; sequence-specific distal enhancer binding RNA polymerase II transcription factor activity; transcription factor binding;