



Anti-EXOG monoclonal antibody, clone 3G8 (DMABT-H16281)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-EXOG Monoclonal Antibody
Target	EXOG
Immunogen	ENDOGL1 partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3G8
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA
Size	100 µg
Preservative	None
Storage	Aliquot and store at -20°C or lower. Avoid repeated freezing and thawing. Multiple freeze/thaw cycles may result in decreased performance.

GENE INFORMATION

Gene Name [EXOG endo/exonuclease \(5-3\), endonuclease G-like \[Homo sapiens \]](#)

Official Symbol	EXOG
Synonyms	EXOG; endo/exonuclease (5-3), endonuclease G-like; ENDOGL1, ENDOGL2, endonuclease G like2, endonuclease G like 2; nuclease EXOG, mitochondrial; ENGL; ENGL a; ENGL b; endo G-like 1; endonuclease G-like 1; endonuclease G-like 2; ENGLA; ENGLB; ENGL-a; ENGL-b; ENDOGL1; ENDOGL2; MGC125944; MGC125945;
Entrez Gene ID	9941
Protein Refseq	NP_001138936
UniProt ID	Q9Y2C4
Chromosome Location	3p21.3
Function	endonuclease activity; hydrolase activity; metal ion binding; nucleic acid binding;