



Anti-Amitriptyline monoclonal antibody, clone 304 (DMABT-Z59177)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-Amitriptyline Monoclonal Antibody
Target	Amitriptyline
Immunogen	Amitriptyline conjugated to KLH.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	304
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Amitriptyline is in a class of drugs called tricyclic antidepressants. Amitriptyline is used to relieve symptoms of depression and also sometimes to treat certain types of pain.
Keywords	AMITRIPTYLINE; AMITRYPTYLLINE; NA; Adepress; Adepril; amitriptilina; Amitriptylin; Amitryptiline; Amitryptyline; Amytriptiline; Damilan; Damilen; Damitriptyline; elani; Elanil; Elavil; Flavyl; Lantron; Laroxil; Laroxyl; Lentizol; mk230; MK-230; N 750