



Anti-Gentamycin monoclonal antibody, clone 204 (DMABT-Z60417)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-Gentamycin Monoclonal Antibody
Target	Gentamycin
Immunogen	Full length protein conjugated to KLH.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	N/A
Clone	204
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.

Keywords

Gentamicin;
