



Anti-AUB (internal region) polyclonal antibody (DPAB-DC1673)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal region of drosophila aub.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Conjugate	Unconjugated
Applications	WB, IP, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	aub aubergine [Drosophila melanogaster (fruit fly)]
Official Symbol	AUB
Synonyms	AUB; aubergine; Aub; AUB; CG6137; Dmel\CG6137; ms(2)32D; sti; sting; CG6137-PA;

CG6137-PC; aub-PA; aub-PC; stellate-interacting gene;

Entrez Gene ID	34524
Protein Refseq	NP_476734
UniProt ID	O76922
Chromosome Location	32C1-32C1; 2-39 cM
Pathway	Dorso-ventral axis formation;
Function	RNA binding; piRNA binding;