



Anti-FASTKD5 (aa 671-764) polyclonal antibody (DPAB-DC2709)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FASTKD5 (FAST kinase domains 5) is a protein-coding gene. Diseases associated with FASTKD5 include schizophrenia. GO annotations related to this gene include protein kinase activity.
Immunogen	FLJ13149 (NP_068598, 671 a.a. ~ 764 a.a) partial recombinant protein with GST tag. The sequence is NVADKSGAMEMAGLCPAACMQTPRMKLAVQFTNRNQYCYGSRDLLGLHNMKRRQLARLGY RVVELSYWEWLPLLKRTREKLAFLEKVFVTSAL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [FASTKD5 FAST kinase domains 5 \[Homo sapiens \(human\) \]](#)

Official Symbol	FASTKD5
Synonyms	FASTKD5; FAST kinase domains 5; dJ1187M17.5; FAST kinase domain-containing protein 5;
Entrez Gene ID	60493
Protein Refseq	NP_068598
UniProt ID	Q7L8L6
Chromosome Location	20p13
Function	poly(A) RNA binding; protein kinase activity;