



Anti-THNSL1 (aa 151-250) polyclonal antibody (DPAB-DC3295)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	THNSL1 (threonine synthase-like 1 (<i>S. cerevisiae</i>)) is a protein-coding gene. Diseases associated with THNSL1 include pneumonia, and alzheimer's disease. GO annotations related to this gene include molecular_function and threonine synthase activity. An important paralog of this gene is THNSL2.
Immunogen	THNSL1 (NP_079114, 151 a.a. ~ 250 a.a) partial recombinant protein with GST tag. The sequence is IVYLDVPLLDLICRLKLMKTDRIVGQNSGTSMDLLKFRRQYYKKWYDARVFCESGASPE EVADKVLNAIKRYQDQVDFSTRHVWPEACEQKVSASF
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [THNSL1 threonine synthase-like 1 \(*S. cerevisiae*\) \[*Homo sapiens* \(human\) \]](#)

Official Symbol	THNSL1
Synonyms	THNSL1; threonine synthase-like 1 (S. cerevisiae); TSH1; threonine synthase-like 1;
Entrez Gene ID	79896
Protein Refseq	NP_079114
UniProt ID	Q8IYQ7
Chromosome Location	10p12.1
Function	molecular_function;