



Anti-C Reactive Protein polyclonal antibody (DPAB-DC4139)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	By immunoelectrophoresis and gel diffusion, This antibody versus fluids with high levels of CRP shows a single band.
Target	C Reactive Protein
Immunogen	Dog CRP.
Source/Host	Goat
Species Reactivity	Dog
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	250 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CRP

GENE INFORMATION

Entrez Gene ID [1401](#)

UniProt ID [P02741](#)
