



## Anti-C Reactive Protein polyclonal antibody (DPAB-DC4139)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	By immunoelectrophoresis and gel diffusion, This antibody versus fluids with high levels of CRP shows a single band.
<b>Target</b>	C Reactive Protein
<b>Immunogen</b>	Dog CRP.
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Dog
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	250 µl
<b>Buffer</b>	In buffer containing 0.1% sodium azide
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

### BACKGROUND

**Keywords** CRP

# GENE INFORMATION

Entrez Gene ID [1401](#)

---

UniProt ID [P02741](#)

---